

Rubisco (Large Chain) Monoclonal Antibody



MAC00378

Product Information

Size:

50ug

Reactivity:

Plants

Source:

Mouse

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:5000

Protein Background:

Gene ID:

Rubisco (Large Chain)

Uniprot

Synonyms:


Immunogen:

Protein

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western blot analysis of Arabidopsis with Rubisco(Large Chain) Mouse mAb diluted at 1)1:2,000 2)1:5,000