

S1-Tag Monoclonal Antibody



MACO0385

Product Information

Size:

50ug

Reactivity:

Source:

Mouse

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

S1-Tag

Uniprot

Synonyms:

Immunogen:

Synthetic Peptide

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western blot analysis of Recombinant S1-Tag Protein with S1-Tag Mouse mAb diluted at 1:5,000.