

TUBE1 Monoclonal Antibody



MACO0397

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Mouse

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:200

Protein Background:**Gene ID:**

TUBE1

Uniprot

Q9UJT0

Synonyms:

Tubulin epsilon chain (Epsilon-tubulin)


Immunogen:

Synthetic Peptide of Epsilon Tubulin


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Immunohistochemical analysis of paraffin-embedded Human Colon Carcinoma Tissue using Epsilon Tubulin Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Colon Carcinoma Tissue using Epsilon Tubulin Mouse mAb diluted at 1:200.