

Glu-Glu Tag Monoclonal Antibody



MAC00426

Product Information

Size:

50ug

Reactivity:

Source:

Mouse

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:5000

Protein Background:

Gene ID:

Glu-Glu Tag

Uniprot

N/A

Synonyms:

Immunogen:

Synthetic Peptide of Glu-Glu Tag

Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images

N/A

N/A