

# AMBP Monoclonal Antibody



MACO0546

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Mouse

**Isotype:**

IgG3

**Applications:**

ELISA, IHC

**Recommended dilutions:**

IHC:1:50-1:500

**Protein Background:**

Inter-alpha-trypsin inhibitor inhibits trypsin, plasmin, and lysosomal granulocytic elastase. Inhibits calcium oxalate crystallization. Trypstatin is a trypsin inhibitor

**Gene ID:**

AMBP

**Uniprot**

P02760

**Synonyms:**

Protein HC, Alpha-1 microglycoprotein, ITI-LC, Bikunin, EDC1, HI-30

**Immunogen:**


Recombinant Human protein AMBP

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemical of paraffin-embedded human liver tissue using MACO0546 at dilution of 1:200