

# Sdc1 Monoclonal Antibody



MACO0556

---

## Product Information

**Size:**

50ug

**Reactivity:**

Rat

**Source:**

Mouse

**Isotype:**

IgG1

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Cell surface proteoglycan that bears both heparan sulfate and chondroitin sulfate and that links the cytoskeleton to the interstitial matrix.

**Gene ID:**

Sdc1

**Uniprot**

P26260

**Synonyms:**

Syndecan-1, SYND1, CD\_antigen: CD138

**Immunogen:**

Recombinant Rat Syndecan-1 protein

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A