

# Dehydroepiandrosterone(DHEA) Monoclonal Antibody

## MACO0588



---

### Product Information

**Size:**

50ug

**Reactivity:****Source:**

Mouse

**Isotype:**

IgG1

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Dehydroepiandrosterone is an important endogenous steroid hormone. It is the most abundant circulating steroid hormone in humans, in whom it is produced in the adrenal glands, the gonads, and the brain, where it functions predominantly as a metabolic intermediate in the biosynthesis of the androgen and estrogen sex steroids.

**Gene ID:**

Dehydroepiandrosterone(DHEA)

**Uniprot****Synonyms:**

DHEA, androstenolone, dehydroisoandrosterone, prasterone, diandron, psicosterone, trans-dehydroandrosterone, diandron, 17-hormoforin

**Immunogen:**

Dehydroepiandrosterone-BSA conjugate

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A