

ACR Antibody



PACO00001

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

ACR

Uniprot

P10323

Synonyms:

ACR; ACRS; Acrosin


Immunogen:

Synthesized peptide derived from the N-terminal region of human Acrosin.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HepG2 cells using Acrosin Polyclonal Antibody.