

Cleaved-CTSV (L114) Antibody



PACO00024

Product Information

Size:

50ug

Reactivity:

Human, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

CTSV

Uniprot

O60911

Synonyms:

CTSL2; CATL2; CTSU; CTSV; Cathepsin L2; Cathepsin U; Cathepsin V


Immunogen:

Synthesized peptide derived from the Internal region of human Cathepsin L2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HeLa cells using Cleaved-Cathepsin L2 (L114) Polyclonal Antibody.