

Cleaved-F2 (R327) Antibody



PACO00040

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

F2

Uniprot

P00734

Synonyms:

F2; Prothrombin; Coagulation factor II


Immunogen:

Synthesized peptide derived from the Internal region of human Thrombin APII.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HeLa cells using Cleaved-Thrombin APII (R327) Polyclonal Antibody.