

Cleaved-MMP1 (F100) Antibody



PACO00054

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

MMP1

Uniprot

P03956

Synonyms:

MMP1; CLG; Interstitial collagenase; Fibroblast collagenase; Matrix metalloproteinase-1; MMP-1


Immunogen:

Synthesized peptide derived from the Internal region of human MMP-1 22k.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of A549 cells using Cleaved-MMP-1 22k (F100) Polyclonal Antibody.