

GAPDH Antibody



PACO00084

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat, Rabbit, Chicken,
Monkey, Sheep, X

Source:

Rabbit

Isotype:

IgG

Applications:

WB, IHC

Recommended dilutions:

IHC:1:200

Protein Background:**Gene ID:**

GAPDH

Uniprot

P04406

Synonyms:

GAPDH; GAPD; CDABP0047; OK/SW-cl.12; Glyceraldehyde-3-phosphate
dehydrogenase; GAPDH; Peptidyl-cysteine S-nitrosylase GAPDH


Immunogen:

Recombinant Protein.

Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images



Western blot analysis of 293T (1), Rat brain (2), NIH 3T3 (3), Sheep Muscle (4), Rabbit testis (5), diluted at 1:20000.