

EFHD1 Antibody



PACO00089

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

WB

Recommended dilutions:

Protein Background:

Gene ID:

EFHD1

Uniprot

Q9BUP0

Synonyms:

EF-hand domain-containing protein D1


Immunogen:

Recombinant Protein.

Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images



Western blot analysis of 293T cell lysate, diluted at 1:2000.