

# Alliinase Antibody



PACO00093

---

## Product Information

**Size:**

50ug

**Reactivity:**

Allium sativum L

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

**Uniprot**

**Synonyms:**

**Immunogen:**


Purified Protein.

**Storage:**

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

## Product Images

---



Western blot analysis of purified alliinase, diluted at 1:5000.