

Phospho-VASP (S157) Antibody



PACO00273

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

VASP

Uniprot

P50552

Synonyms:

VASP; Vasodilator-stimulated phosphoprotein; VASP


Immunogen:

Synthesized peptide derived from human VASP around the phosphorylation site of S157.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COS7 cells using Phospho-VASP (S157) Polyclonal Antibody.