

ACTBL2 Antibody



PACO00366

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

ACTBL2

Uniprot

Q562R1

Synonyms:

ACTBL2; Beta-actin-like protein 2; Kappa-actin


Immunogen:

Synthesized peptide derived from the Internal region of human Actin-Kappa.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of A549 cells using Actin-Kappa Polyclonal Antibody.