

ADCK4 Antibody



PACO00371

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:**Gene ID:**

COQ8B

Uniprot

Q96D53

Synonyms:

ADCK4; Uncharacterized aarF domain-containing protein kinase 4


Immunogen:

Synthesized peptide derived from the N-terminal region of human ADCK4.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HepG2 cells using ADCK4 Polyclonal Antibody.