

# AIG1 Antibody



PACO00379

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:****Gene ID:**

AIG1

**Uniprot**

Q9NVV5

**Synonyms:**

AIG1; CGI-103; Androgen-induced gene 1 protein; AIG-1

**Immunogen:**


Synthesized peptide derived from the C-terminal region of human AIG-1.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of Jurkat cells using AIG-1 Polyclonal Antibody.