

ANXA10 Antibody



PACO00396

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

ANXA10

Uniprot

Q9UJ72

Synonyms:

ANXA10; ANX14; Annexin A10; Annexin-10; Annexin-14

Immunogen:

Synthesized peptide derived from the C-terminal region of human Annexin X.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of NIH-3T3 cells using Annexin X Polyclonal Antibody.

