

TFAP2C Antibody



PACO00402

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

TFAP2C

Uniprot

Q92754

Synonyms:

TFAP2C; Transcription factor AP-2 gamma; AP2-gamma; Activating enhancer-binding protein 2 gamma; Transcription factor ERF-1

Immunogen:

Synthesized peptide derived from the C-terminal region of human AP-2gamma.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of HuvEc cells using AP-2gamma Polyclonal Antibody.

