

# APOL2 Antibody



PACO00409

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300,  
IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

APOL2

**Uniprot**

Q9BQE5

**Synonyms:**

APOL2; Apolipoprotein L2; Apolipoprotein L-II; ApoL-II

**Immunogen:**


Synthesized peptide derived from the Internal region of human ApoL2.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of COLO205 cells using ApoL2 Polyclonal Antibody.