

# MIP Antibody



PACO00411

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300

**Protein Background:****Gene ID:**

MIP

**Uniprot**

P30301

**Synonyms:**

MIP; AQP0; Lens fiber major intrinsic protein; Aquaporin-0; MIP26; MP26

**Immunogen:**

Synthesized peptide derived from the Internal region of human AQP0.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

Western Blot analysis of HT29 cells using AQP0 Polyclonal Antibody.

