

ARMC6 Antibody



PACO00426

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

ARMC6

Uniprot

Q6NXE6

Synonyms:

ARMC6; Armadillo repeat-containing protein 6

Immunogen:

Synthesized peptide derived from the C-terminal region of human ARMC6.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of HT29 cells using ARMC6 Polyclonal Antibody.

