

ARRB1 Antibody



PACO00429

Product Information

Size:

50ug

Reactivity:

Human, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:

Gene ID:

ARRB1

Uniprot

P49407

Synonyms:

ARRB1; ARR1; Beta-arrestin-1; Arrestin beta-1


Immunogen:

Synthesized peptide derived from human Arrestin-Beta-1 around the non-phosphorylation site of S412.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of A549 cells using Arrestin-Beta-1 Polyclonal Antibody.