

ATN1 Antibody



PACO00445

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

ATN1

Uniprot

P54259

Synonyms:

ATN1; D12S755E; DRPLA; Atrophin-1; Dentatorubral-pallidolusian atrophy protein

Immunogen:

Synthesized peptide derived from the N-terminal region of human Atrophin-1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of NIH-3T3 cells using Atrophin-1 Polyclonal Antibody.

