

BLK Antibody



PACO00468

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:

Gene ID:

BLK

Uniprot

P51451

Synonyms:

BLK; Tyrosine-protein kinase Blk; B lymphocyte kinase; p55-Blk


Immunogen:

Synthesized peptide derived from human Blk around the non-phosphorylation site of Y501.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of 293T cells using Blk Polyclonal Antibody.