

C1QL2 Antibody



PACO00486

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

C1QL2

Uniprot

Q7Z5L3

Synonyms:

C1QL2; Complement C1q-like protein 2


Immunogen:

Synthesized peptide derived from the C-terminal region of human C1qL2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HuvEc cells using C1qL2 Polyclonal Antibody.