

CA14 Antibody



PACO00492

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

CA14

Uniprot

Q9ULX7

Synonyms:

CA14; Carbonic anhydrase 14; Carbonate dehydratase XIV; Carbonic anhydrase XIV; CA-XIV

Immunogen:

Synthesized peptide derived from the Internal region of human CA XIV.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of 293 cells using CA XIV Polyclonal Antibody.

