

CCDC102B Antibody



PACO00538

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

CCDC102B

Uniprot

Q68D86

Synonyms:

CCDC102B; C18orf14; Coiled-coil domain-containing protein 102B


Immunogen:

Synthesized peptide derived from the Internal region of human CCDC102B.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COLO205 cells using CCDC102B Polyclonal Antibody.