

# CLDN6 Antibody



PACO00608

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300,  
IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

CLDN6

**Uniprot**

P56747

**Synonyms:**

CLDN6; Claudin-6; Skullin

**Immunogen:**


Synthesized peptide derived from the Internal region of human Claudin-6.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of Jurkat cells using Claudin-6 Polyclonal Antibody.