

CYP1A1/CYP1A2 Antibody



PACO00661

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:**Gene ID:**

CYP1A1/CYP1A2

Uniprot

P04798/P05177

Synonyms:

CYP1A1; Cytochrome P450 1A1; CYP1A1; Cytochrome P450 form 6; Cytochrome P450-C;
Cytochrome P450-P1; CYP1A2; Cytochrome P450 1A2; CYP1A2; Cytochrome P(3)450;
Cytochrome P450 4; Cytochrome P450-P3


Immunogen:

Synthesized peptide derived from the Internal region of human CYP1A1/2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HT29 cells using CYP1A1/2 Polyclonal Antibody.