

CYP2A13 Antibody



PACO00663

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:

Gene ID:

CYP2A13

Uniprot

Q16696

Synonyms:

CYP2A13; Cytochrome P450 2A13; CYP11A13


Immunogen:

Synthesized peptide derived from the Internal region of human CYP2A13.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of A549 cells using CYP2A13 Polyclonal Antibody.