

CYP2S1 Antibody



PACO00667

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:

Gene ID:

CYP2S1

Uniprot

Q96SQ9

Synonyms:

CYP2S1; Cytochrome P450 2S1; CYP11S1


Immunogen:

Synthesized peptide derived from the Internal region of human CYP2S1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COLO205 cells using CYP2S1 Polyclonal Antibody.