

CYP4Z2P Antibody




PACO00672

Product Information

Size: 50ug	Protein Background:
Reactivity: Human	Gene ID: CYP4Z2P
Source: Rabbit	Uniprot Q5VVE4/Q8N1L4
Isotype: IgG	Synonyms: Cytochrome P450 4Z1/Putative cytochrome P450 family member 4Z2
Applications: ELISA, WB, IHC	Immunogen: Synthesized peptide derived from the Internal region of human CYP4Z1/2.
Recommended dilutions: IHC:1:100-1:300	Storage: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of Jurkat RAT-MUSCLE cells using CYP4Z1/2 Polyclonal Antibody.