

DDX51 Antibody



PACO00693

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

DDX51

Uniprot

Q8N8A6

Synonyms:

DDX51; ATP-dependent RNA helicase DDX51; DEAD box protein 51


Immunogen:

Synthesized peptide derived from the C-terminal region of human DDX51.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HepG2 cells using DDX51 Polyclonal Antibody.