

POLA1 Antibody



PACO00714

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

POLA1

Uniprot

P09884

Synonyms:

POLA1; POLA; DNA polymerase alpha catalytic subunit; DNA polymerase alpha catalytic subunit p180


Immunogen:

Synthesized peptide derived from the N-terminal region of human DNA pol alpha.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of Jurkat cells using DNA POLL Ipha Polyclonal Antibody.