

EEF1G Antibody



PACO00739

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

EEF1G

Uniprot

P26641

Synonyms:

EEF1G; EF1G; Elongation factor 1-gamma; EF-1-gamma; eEF-1B gamma


Immunogen:

Synthesized peptide derived from the Internal region of human EF-1 gamma.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HepG2 cells using EF-1 gamma Polyclonal Antibody.