

# ZC3H8 Antibody



PACO00808

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

ZC3H8

**Uniprot**

Q8N5P1

**Synonyms:**

ZC3H8; ZC3HDC8; Zinc finger CCCH domain-containing protein 8

**Immunogen:**


Synthesized peptide derived from the Internal region of human Fliz1.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of HepG2 cells using Fliz1 Polyclonal Antibody.