

FGFR1OP Antibody



PACO00814

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:

Gene ID:

FGFR1OP

Uniprot

O95684

Synonyms:

FGFR1OP; FOP; FGFR1 oncogene partner


Immunogen:

Synthesized peptide derived from the C-terminal region of human FOP.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HepG2 cells using FOP Polyclonal Antibody.