

FOXD3 Antibody



PACO00818

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

FOXD3

Uniprot

Q9UJU5

Synonyms:

FOXD3; HFH2; Forkhead box protein D3; HNF3/FH transcription factor genesis


Immunogen:

Synthesized peptide derived from the Internal region of human FoxD3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of 293 cells using FoxD3 Polyclonal Antibody.