

FZD9 Antibody



PACO00825

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

FZD9

Uniprot

O00144

Synonyms:

FZD9; FZD3; Frizzled-9; Fz-9; hFz9; FzE6; CD antigen CD349


Immunogen:

Synthesized peptide derived from the C-terminal region of human Frizzled-9.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COS7 cells using Frizzled-9 Polyclonal Antibody.