

FUK Antibody



PACO00827

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:**Gene ID:**

FUK

Uniprot

Q8N0W3

Synonyms:

FUK; L-fucose kinase; Fucokinase


Immunogen:

Synthesized peptide derived from the N-terminal region of human Fucokinase.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COS7 cells using Fucokinase Polyclonal Antibody.