

GPR116 Antibody



PACO00865

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:

Gene ID:

ADGRF5

Uniprot

Q8IZF2

Synonyms:

GPR116; KIAA0758; Probable G-protein coupled receptor 116


Immunogen:

Synthesized peptide derived from the N-terminal region of human GPR116.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of A549 cells using GPR116 Polyclonal Antibody.