

GPR50 Antibody



PACO00877

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:**Gene ID:**

GPR50

Uniprot

Q13585

Synonyms:

GPR50; Melatonin-related receptor; G protein-coupled receptor 50; H9


Immunogen:

Synthesized peptide derived from the Internal region of human GPR50.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of A549 cells using GPR50 Polyclonal Antibody.