

GRB14 Antibody



PACO00884

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

GRB14

Uniprot

Q14449

Synonyms:

GRB14; Growth factor receptor-bound protein 14; GRB14 adapter protein


Immunogen:

Synthesized peptide derived from the Internal region of human GRB14.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of MCF7 cells using GRB14 Polyclonal Antibody.