

GRTP1 Antibody



PACO00890

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

GRTP1

Uniprot

Q5TC63

Synonyms:

GRTP1; TBC1D6; Growth hormone-regulated TBC protein 1; TBC1 domain family member 6


Immunogen:

Synthesized peptide derived from the Internal region of human GRTP1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of 293 cells using GRTP1 Polyclonal Antibody.