

HOXD12 Antibody



PACO00933

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

HOXD12

Uniprot

P35452

Synonyms:

HOXD12; HOX4H; Homeobox protein Hox-D12; Homeobox protein Hox-4H


Immunogen:

Synthesized peptide derived from the C-terminal region of human HoxD12.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of K562 cells using HoxD12 Polyclonal Antibody.