

IL21 Antibody



PACO00960

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

IL21

Uniprot

Q9HBE4

Synonyms:

IL21; Interleukin-21; IL-21; Za11


Immunogen:

Synthesized peptide derived from the Internal region of human IL-21.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of Jurkat cells using IL-21 Polyclonal Antibody.