

INHBC Antibody



PACO00967

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

INHBC

Uniprot

P55103

Synonyms:

INHBC; Inhibin beta C chain; Activin beta-C chain


Immunogen:

Synthesized peptide derived from the Internal region of human Inhibin Beta-C.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of A549 cells using Inhibin Beta-C Polyclonal Antibody.