

IQGAP3 Antibody



PACO00980

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

IQGAP3

Uniprot

Q86VI3

Synonyms:

IQGAP3; Ras GTPase-activating-like protein IQGAP3


Immunogen:

Synthesized peptide derived from the N-terminal region of human IQGAP3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of Jurkat cells using IQGAP3 Polyclonal Antibody.